

# ENDOSCOPE VIDEO CAMERA



## Characteristics

- Color management
- Automatic, manual, preset white balance setting
- Control of the operating mode of the image quality improvement system
- Full HD video recording on SD card-optional zoom
- Adapter with variable focal length- optional.



## Technical specification

Maximum Resolution:	1920 x 1080
Aspect Ratio:	16:9
Sensitivity:	Min 0.1 lx
Signal to Noise Ratio:	>48 dB
Video Output:	DVI/HDMI (DVI-D) 3G-SDI (BNC) CVC (BNC-Composite)