




VIDEO OTOSCOPE





Sensor with precise optic



Technical Specification

Diameter Ø (Probe/Ear Funnel):	4.0mm/2.5mm + 4.0mm
Focal Distance:	8-20 mm
Interfaces	USB
Power supply/consumption:	USB/ max. 0.5 A
Risc Class (according to medical law):	1
Protection Class (against electric shock):	class II
Operating Mode:	Continuous operation
Protection Class:	IP 65
Temperature (Transport/operation):	-10° C to +70° C/ 0° C to 40° C
Relative Humidity:	0% to 95 %
Air Pressure:	950hpa to 1050hpa

